

AUTOMATED, LARGE SCALE TARGETED DIFFERENTIATION OF iPS CELLS FOR CELL THERAPY APPLICATIONS

Itsuki Kinoshita¹, Yohei Hayashi², Yoshiki Nakashima², Takaki Yamaguchi¹, Nathan D Frank³, Mindy Miller³, Masatsugu Igarashi¹, Masayoshi Tsukahara²

1: Terumo Corporation, 2: CiRA Foundation, 3: Terumo BCT, Inc.



Introduction

Induced pluripotent stem (iPS) cells offer great potential as an unlimited source for regenerative therapies. However, their clinical use remains constrained by two key challenges: achieving efficient differentiation into target cell types and scaling production to clinically relevant quantities. Here, we present a large-scale protocol for generating neural stem cells (NSC) from iPS cells using the Quantum Flex system (Quantum, Terumo BCT).

Materials and Methods

[Figure 1: Culture Process]

Day -1: Coating

Day 0: Load/Attach cells

Day 1: Start Differentiation

Day 1-10: Differentiation


Day 10: Harvest

- Coating Reagent: Laminin 521 (Bioamina, 1 mg)
- Seed Density: 1300 cells/cm² (total 2.6 x 10⁶ cells)
- Medium: mTeSR Plus (STEMCELL Technologies)
- Medium Change into Differentiation Medium
Differentiation Media Composition:
StemFit AK03N (Ajinomoto, without C solution)
+ 2 μM SB431542 + 2 μM DMH1 (FUJIFILM Wako, 037-24293, 041-33881)
- Perform Medium Change every other day using new protocol (Figure 2)
- Automated Metabolism Monitoring by Bioprofile FLEX2 (FLEX2, Nova Biomedical) integrated with Quantum Flex
- Dissociation Reagent: Accutase (60 mL, 4 min)

[Figure 2: Medium Replacement Protocol]

Step 1:

Extracapillary circuit



11 mL/min

Step 2:

Intracapillary circuit (reactor exterior)

10 mL/min



Step 3:

Intracapillary circuit (reactor interior)



0.3 mL/min


This new medium-exchange process achieved:

- ✓ Minimal cell loss during medium exchange (<1.5%)
- ✓ High medium-exchange efficiency (>90%)
- ✓ Rapid medium exchange (3.5 h)


Results

[Figure 3: Monitoring by FLEX2]

Glucose and Lactate readings




pH



- ✓ Integration with FLEX2 enabled high-frequency auto monitoring that **simplified culture condition management**


[Figure 4: mRNA Expression]



(n=3, t-test, *: P < 0.05, **: P < 0.005, ***: P < 0.0005, NS = not significant)

- ✓ Confirmed **loss of pluripotency markers and induction of neural stem cell markers**
- ✓ Achieved **more consistent differentiation efficiency with the Quantum** compared with manual culture

[Figure 5: Flow cytometry]



Conclusion

We established a robust and scalable protocol for differentiating iPS cells into neural stem cells using the Quantum Flex system. This platform lays the groundwork for expanding future differentiation applications and advancing clinically relevant cell manufacturing.