



◆ QHJ14s03 (iPS cells expressing the highest HLA in Japan※1)

Clone ID	QHJ14s03	Product	Human iPS cells
Source	Peripheral Blood, Human	Race	Japanese
Passage No.	10	Gender	Male
Lot No.	20150802-02	Manufacture Dates	Aug. 2 nd , 2015
Culture medium	StemFit AK03N	Substrate	iMatrix-511MG
Culture Method	Feeder-free (※2)		
Plasmids for reprogramming	pCE-hSK, pCE-hUL, pCE-hOCT3/4, pCE-mp53DD, pCXB-EBNA1		
Use and Provision of this cell stock	Please check our web site ; https://www.cira-foundation.or.jp/e/project/stock.html		

(※1) Reference; Okita, *et. al.*, Nat Methods. 2011 8(5): 409-412

(※2) Reference; Nakagawa, *et. al.*, Nat Biotechnol. 2008 26(1):101-106

Test Result

Test	Method	Result
Sterility	BacT/ALERT	Negative
Mycoplasma	PCR	Negative
Endotoxin	LAL	≤ 5 EU/mL
Virus (HBV, HCV, HIV, HTLV, Parvovirus B19)	PCR	Negative
HLA typing (HLA-A, B, DR)	PCR-SBT	Consistent with the donor cells
STR genotyping	PCR	Consistent with the donor cells
Morphology	Microscope	Consistent with human ES cells
Karyotype	Conventional Giemsa analysis G-banding	46,XY[20]
Plasmid remnants	qPCR	Below the limit of quantification
CNV(※3)	WGS, SNP	No de novo CNVs (>1kbp) were found in COSMIC Cancer Gene Census (ver.74) and Shibata list(※5) .
SNV/Indel(※4)	WGS, WES	No de novo non-synonymous SNVs/Indels were found in COSMIC Cancer Gene Census (ver.74) and Shibata list(※5) .
Undifferentiated markers	Microarray(※7)	<i>POU5F1</i> : 5.1%、 <i>NANOG</i> : 9.3 % (Relative expression levels of <i>GAPDH</i>)
	Flow cytometry (※7)	TRA-1-60: 94.2% SSEA4: 99.9% TRA-2-49: 100.0 %



Thawed postnatal cells	Counting the number of the cells (※6, 7)	1.25×10^5 cells (Survival rate; 94%)
Number of proliferating cells after thawing	Counting the number of the cells after culturing for 9 days(※6, 7).	15.4×10^5 cells (Number of seeded cells; 6.5×10^4 cells)
Doubling time (h)	Counting the number of the cells (※6, 7)	P13→P14: 36.8 P14→P15: 23.2 P15→P16: 29.1 P16→P17: 25.7 P17→P18: 24.9

(※3) CNV; Copy Number Variation


(※4) SNV/Indel; Single nucleotide variants /Insertion Deletion

(※5) The PMDA Science Board “Current Perspective on Evaluation of Tumorigenicity of Cellular- and Tissue-based Products Derived from induced Pluripotent Stem Cells (iPSCs) and iPSCs as Their Starting Materials” (Cellular- and Tissue-based Products Subcommittee, 20 August 2013)

(※6) ThermoFisher Countess®

(※7) The result of # 1 out of 3 frozen stocks is shown.

■ Image



Please contact us if you have any questions.

ips-request@cira-foundation.or.jp



Reprint or reproduction of this page without permission is prohibited.